



ENQUIRY REF. No.	KAPL/ME/013/1066
DATE	12/08/2024
DUE DATE	14/08/2024 TILL: 13.00 HRS

Dear Sir,

Please submit your lowest and competitive offer in a SEALED ENVELOPE, DULY SUPERSCRIBING OUR ABOVE ENQUIRY REF. NO., DATE and DUE DATE on it/ OR MAIL, with other details of F.O.R terms, Taxes, Credit period, Delivery offered, Name of the Make, Detailed Specification etc., for below mentioned material/s

SL. NO.	ITEM CODE	ITEM DESCRIPTION	UOM	QTY
1	MIV3001	VHP SPRAY VALVE FOR ISOLATOR MDL:SMPL-S150 SNOWBE	NOS	01

Please ensure that your offer reaches us on or before Due Date by courier OR Speed post OR you can also mail us to our email: puren@kaplindia.com

OTHER TERMS:

- | | |
|---------------------------------|------------------|
| 1. F.O.R TERMS | : DOOR DELIVERY |
| 2. GST % | : PLEASE SPECIFY |
| 3. PACKING & FORWARDING CHARGES | : NOT APPLICABLE |
| 4. CREDIT PERIOD | : 30 DAYS |
| 5. DELIVERY OFFERED | : PLEASE SPECIFY |

NOTE: IN CASE YOU ARE NOT QUOTING PLEASE SEND THE REGRET LETTER.

Thanking you,

Yours faithfully,
For KARNATAKA ANTIBIOTICS
& PHARMACEUTICALS LIMITED


YUVARAJA M
DEPUTY MANAGER PURCHASE DEPT

ENGINEERING AND PROJECT DEPARTMENT

MATERIAL REQUISITION

PR# 90241-02565

Date : 05/03/2024

Very Urgently Required.

stores

Through : HOD - E&P

Sl. No. 688

Machine Name and Details : Isolator

Location : OPP3 & Microbiology

Following are the details :

Sl. No.	Item Code No.	Item Description	UOM	Qty.	Required Date	P.R. No. & Date	Remarks
1.	MIV 3001	VHP Spray valve for ISOLATOR Model: SMPL-S150 make: Snow Bell machines Pvt. Ltd.	no.	01	25/03/2024		Required for OPP3 & C Microbiology Sampling Isolator. * Service / MOM report attached for reference
2.							

M/BA vendor

6) Door Lock Toggle clamps were loosened; Tightening is done.

7) VHP Spray Valve issue Checked and found faulty: Tube is blocked in order to arrest the leakage.

Required Components for replacement:-

- 1) Door Gasket (10mtrs).
- 2) Pneumatic butter fly valve (1no).
- 3) Festo Solenoid Valve (1no).
- 4) VHP Spray Valve (1no).
- 5) LED Light strip (6mtrs).
- 6) Toggle clamps(6nos).
- 7) Pneumatic Door Cylinder (1no).

Sampling Isolator Chamber's all leakages are arrested. The Pressure Hold & Leak Test is passing in Auto mode and reports are printed in presence of E&P and QC Microbiology departments Executives. The unit is under observation. KAPL will indent spares. After receiving the spares, the servicing will be carried out.